

CHARIOT®

Protected by
U.S. Patent No. D668,422
Patent Pending



It's a Lift on Wheels!



***Transport your scooter or
powerchair safely with a
full compliment of cars,
SUV's, CUV's and Mini Vans.
Incredible!***

 **BRUNO®**
for your independence

It's a Lift on Wheels!

A Most Remarkable Development!

Many people who depend on scooters and power wheelchairs for mobility have reluctantly traded in their small or mid-size cars for larger and higher vehicles to transport their personal mobility devices.

Now there's an alternative! Introducing the Chariot — an incredible new “lift on wheels” that empowers a small vehicle to handle a load of up to 350 lb with ease.

Positioned behind a vehicle equipped with a frame mounted Class I or II receiver hitch, the Chariot uses two independent wheels that swivel — following the vehicle precisely. With its own suspension, this remarkable lift virtually eliminates any additional weight or force on the vehicle!

It's a miracle!

Look At What People Are Saying About the Chariot...

“It's great! I don't feel the vehicle bounce like with other common lifts.”

“It's so easy to drive...the lift follows right along, even in reverse.”


“I think it's very nice. You can't even feel it back there.”

“It's not like a trailer..so simple to maneuver!”


“It's just like an extension of my car!”




The Chariot Works On Many Different Cars!




When used on a small car, it's shorter than a Full-Size Pickup Truck



Minivan



Chariot On
Small Car



Full-Size
Pickup Truck

Keep All Your Vehicle Seats!



Chariot pictured with a 2013 Ford Edge featuring Bruno's Valet® Plus passenger seat.




Fits In Most Garages!

NOW AVAILABLE!

for select SUVs, CUVs, Minivans & Wagons.

Note: Must meet the specifications for installing on a hatchback type vehicle

The Chariot Goes Where You Go... Just Drive Normally!



NO JACKKNIFING

Swivel wheels on the Chariot means the lift follows you precisely – No "Jackknifing" – You might not even know it's there.

The Chariot Does Not Change Position!

Chariot Features...



**Self-Tensioning Retractable
Securement Belts -
One Hand Operation**



**Swivel Wheels With
Independent Suspension**



Easy Fold Platform



**Drive On/Off Platform From
Either Side - No Reversing**



Self-Leveling Platform



Optional Corner Markers



**Running/Brake Lights and
Turn Signals**



Manual Back-Up



**Lighted License
Plate Holder
(for use in some states/
provinces)**



**One Button Operation
with On/Off Key**



Optional Spare Tire

**The Chariot Works With
Many Types Of Mobility
Devices...**



Three or Four Wheeled Scooters



Mid-Wheel Drive Powerchairs



**Front-Wheel/Rear-Wheel
Drive Powerchairs**



**Manual Wheelchairs with
Seating Systems**

KEY FEATURES

	CHARIOT
Maximum Weight Capacity	350 lbs/159 kg
Power Up/Down	yes
Manual Assist Fold	yes
Single Switch Operation	yes
Adjustable Platform	yes
Powered By Vehicle Battery	yes
Emergency Backup	yes
Compatible w/Class I or Class II receiver hitch (1 1/4" square) - frame mount only	yes
Uses Flat-4 Wiring Connector (Required)	yes
Optional Spare Tire Available	yes
3 Year Limited Warranty	yes
Can Reinstall on Another Vehicle*	yes
Eligible for Auto Manufacturer Rebate*	yes

***IMPORTANT:** Contact your authorized Bruno dealer for specific details on your vehicle!

See The Chariot In Action...



www.bruno.com



Represented by:

Other great Bruno products...



ELITE STAIRLIFT



CUSTOM CURVED RAIL STAIRLIFT



VERTICAL PLATFORM LIFT



CURB-SIDER® VEHICLE LIFT



OUT-SIDER® VEHICLE LIFT



VALET® SIGNATURE SEATING

BRUNO INDEPENDENT LIVING AIDS, INC.®

1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

