



America's Favorite  
Power Chair



|   | Jazzy®<br>Elite ES-1  | Jazzy®<br>Elite ES Portable   | Jazzy®<br>Elite ES  | Jazzy®<br>Select Elite  | Jazzy®<br>Select 6  | Jazzy®<br>Elite 14  | Jazzy®<br>600 ES  | Jazzy®<br>1113 ATS  | Jazzy®<br>Sport 2   | Jazzy®<br>Elite HD  | Jazzy®<br>614 HD  | Jazzy®<br>1450  |
|---|---|---|---|---|---|---|---|---|---|---|---|---|
| Weight Capacity                           | 300 lbs.  | 300 lbs.  | 300 lbs.  | 300 lbs.  | 300 lbs.  | 300 lbs.  | 300 lbs.  | 300 lbs.  | 300 lbs.  | 450 lbs. <sup>12</sup>  | 450 lbs.  | 600 lbs.  |
| Turning Radius <sup>1</sup>               | 23.5"   | 25.5"   | 24.75"  | 24.75"  | 21.87"  | 24"   | 20.5"   | 20.5"   | 24.75"  | 24"   | 22.75"  | 26.5"   |
| Width <sup>1</sup>                        | 22.75"  | 22.75"  | 22.75"  | 22.75"  | 23.5"   | 24.25"  | 24.5"   | 23"   | 23"   | 24.25"  | 25.125"   | 29.25"  |
| Length <sup>1,2</sup>                     | 34.75"  | 34.75"  | 34.75"  | 34.75"  | 35.375"   | 35.5"   | 35.5"   | 41.5"   | 34.75"  | 35.5"   | 35.75"  | 38.25"  |
| Maximum Speed <sup>3</sup>                | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4 mph   | Up to 4.3 mph   | Up to 4.5 mph   |
| Ground Clearance <sup>1</sup>             | 0.75"   | 1.5"  | 1.5"  | 1.5"  | 1.75"   | 3"  | 2.75"   | 2.5"  | 1.5"  | 3"  | 2.6"  | 1.5"  |
| Front Wheels                              | 5" solid anti-tips  | 3" solid anti-tips  | 3" solid anti-tips  | 3" solid anti-tips  | 6" solid casters  | 3" solid anti-tips  | 6" solid OMNI-casters   | 6" solid anti-tips  | 3" solid casters  | 3" solid anti-tips  | 6" solid casters  | 4" solid anti-tips  |
| Drive Wheels                              | 9" solid  | 9" solid  | 9" solid  | 9" solid  | 10" solid   | 14"   | 14"   | 10" <sup>8</sup>  | 9"  | 14" solid   | 14"   | 14" solid   |
| Rear Wheels                               | 6" solid casters  | 6" solid casters  | 6" solid casters  | 6" solid casters  | 6" solid casters  | 6" solid casters  | 6" solid OMNI-casters   | 8" solid casters  | 6" solid casters  | 6" solid casters  | 6" solid casters  | 9" casters  |
| Drivetrain                                | Two-motor, in-line, front-wheel drive   | Two-motor, in-line, front-wheel drive   | Two-motor, in-line, front-wheel drive   | Two-motor, in-line, front-wheel drive   | Two-motor, in-line, mid-wheel drive   | Two-motor, front-wheel drive  | Two-motor, mid-wheel drive  | Two-motor, mid-wheel drive  | Two-motor, in-line, mid-wheel drive   | Two-motor, front-wheel drive  | Two-motor, Mid-Wheel 6 <sup>®</sup>   | Two-motor, front-wheel  |
| Braking System                            | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   | Regenerative and electro-mechanical   |
| Suspension Type                           | Limited   | Limited   | Limited   | Limited   | Active-Trac <sup>®</sup> (ATX)  | Limited   | Active-Trac <sup>®</sup> (ATX)  | Active-Trac <sup>®</sup> (ATX)  | Limited   | Limited   | Limited   | Swing Arm   |
| High-Back Seat (max. size)                | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 24"W x 24"D   | 24"W x 24"D   | 32"W x 28"D   |
| High-Back Seat-to-Floor Range             | 21.5" - 23.5"   | 21.5" - 23.5"   | 21.5" - 23.5"   | 21.5" - 23.5"   | 21.25" - 23.25"   | 20" - 22"   | 21" - 24"   | 19.25" - 20.75"   | 21.5" - 23.5"   | 20" - 22"   | 20" - 22"   | 20.5" - 22.5"   |
| Specialty Seating (max. size)             | N/A   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | N/A   | N/A   | N/A   |
| Specialty Seat-to-Floor Range             | N/A   | 16.5" - 18.5"   | 16.5" - 18.5"   | 16.5" - 18.5"   | 16.5" - 18.5"   | 16.25" - 18.25"   | 17" - 20"   | 14.75" - 16.5"  | 16.5" - 18.5"   | N/A   | N/A   | N/A   |
| Synergy <sup>®</sup> Seating              | N/A   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 20"W x 20"D   | 24"W x 24"D   | 24"W x 24"D   | 32"W x 28"D   |
| Synergy <sup>®</sup> Seat-to-Floor Range  | N/A   | 17.25" - 19.25"   | 17.25" - 19.25"   | 17.25" - 19.25"   | 16.5" - 18.5"   | 15.75" - 17.75"   | 17.25" - 19.25"   | 15.125" x 17.125"   | 17.25" - 19.25"   | 15.75" - 17.75"   | 15.75" - 17.75"   | 17.25" - 19.25"   |
| Power Elevating Seat Option <sup>11</sup> | N/A   | N/A   | N/A   | N/A   | Yes   | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   |
| TRU-Balance <sup>®</sup>                  | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   | Tilt or Recline <sup>7</sup>  | N/A   | N/A   | N/A   | N/A   |
| Specialty Controls                        | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   | Yes   | N/A   | N/A   | N/A   | Yes   |
| Standard Electronics                      | 40A, PG GC3   | 40A, PG GC3   | 40A, PG GC3   | 40A, PG GC3   | 40A, PG GC3 <sup>10</sup>   | 50A, Dynamic Shark  | 50A, Dynamic Shark  | 50A, PG VSI   | 40A, PG GC3   | 60A, Dynamic Shark  | 50A, Dynamic Shark  | 90A, PG VR2   |
| Battery Size (2 required)                 | 12 volt, U-1  | 12 volt, U-1  | 12 volt, U-1  | 12 volt, U-1  | 12 volt, U-1  | 12 volt, NF-22  | 12 volt, NF-22  | 12 volt, U-1  | 12 volt, U-1  | 12 volt, NF-22  | 12 volt, NF-22  | 12 volt, Group 24   |
| Standard Battery Charger                  | 3A, off-board   | 3A, off-board   | 3A, off-board   | 3A, off-board   | 3.5A, off-board   | 5A, off-board   | 5A, off-board   | 4A, on-board  | 3A, off-board   | 5A, off-board   | 5A, off-board   | 8A, off-board   |
| Per-Charge Range <sup>3,6</sup>           | Up to 12 miles  | Up to 11 miles  | Up to 15 miles  | Up to 15 miles  | Up to 15 miles  | Up to 9.8 miles   | Up to 16.65 miles   | Up to 10.9 miles  | Up to 15 miles  | Up to 9.8 miles   | Up to 11 miles  | Up to 12.14 miles   |
| Battery Weight <sup>9</sup>               | 24.5 lbs. each  | 24.5 lbs. each  | 24 lbs. each  | 24.5 lbs. each  | 24 lbs. each  | 37 lbs. each  | 37 lbs. each  | 24.5 lbs. each  | 24.5 lbs. each  | 37 lbs. each  | 37 lbs.   | 52 lbs. each  |
| Weight of Heaviest Piece <sup>5</sup>     | 68.5 lbs.   | 47.6 lbs.   | 68.5 lbs.   | 68.5 lbs.   | 100 lbs.  | 98.5 lbs.   | 123.6 lbs.  | 77 lbs.   | 68.5 lbs.   | 98.5 lbs.   | 121.5 lbs.  | 195 lbs.  |
| Standard Seat Weight                      | 43 lbs.<br>(comfort, high-back)   | 43 lbs.<br>(comfort, high-back)   | 43 lbs.<br>(comfort, high-back)   | 49 lbs.<br>(comfort, high-back)   | 43 lbs.<br>(comfort, high-back)   | 51.5 lbs.<br>(high-back)  | 51.5 lbs.<br>(high-back)  | 37 lbs.<br>(medium-back)  | 49 lbs.<br>(comfort, high-back)   | 51.5 lbs.<br>(high-back)  | 51.5 lbs.<br>(high-back)  | 51.5 lbs.<br>(high-back)  |
| Standard Seat Size                        | 18" x 18"   | 18" x 18"   | 18" x 18"   | 18" x 18"-20"   | 18" x 18"-20"   | 20" x 20"-18"   | 18" x 18" / 20" x 20"-18"   | 18" x 18"   | 18" x 18"-20"   | 20" x 18"-20"   | 20" x 18"-20"   | 24" x 24"   |
| Medicare Code Coverage                    | <b>Group 1:</b><br>Standard (K0816)   | <b>Group 2:</b> Standard,<br>Portable (K0820/K0821)   | <b>Group 2:</b> Standard<br>(K0822/K0823)   | <b>Group 2:</b> Standard<br>(K0822/K0823)   | <b>Group 2:</b> Standard<br>(K0822/K0823)   | <b>Group 2:</b> Standard<br>(K0822/K0823)   | <b>Group 2:</b> Standard<br>(K0822/K0823)   | <b>Group 2:</b> Single Power<br>(K0835/K0836)   | N/A   | <b>Group 2:</b> Heavy Duty<br>(K0824/K0825)   | <b>Group 2:</b> Heavy Duty<br>(K0824/K0825)   | <b>Group 2:</b> Very Heavy Duty<br>(K0826/K0827), Single Power<br>Very Heavy Duty (K0839)   |
| Warranty                                  | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | Lifetime limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | 5-year limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | Lifetime limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors | <b>K0826/K0827 Models:</b><br>Lifetime limited on frame;<br>13-month limited on<br>electronics; 13-month<br>on drive motors<br><br><b>K0839 Model:</b> Lifetime<br>limited on frame; 2-year<br>limited on electronics;<br>18-month limited on drive<br>motors |

FOOTNOTES:

- <sup>1)</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

<sup>2)</sup> Without front riggings.

<sup>3)</sup> Varies with user weight, terrain, battery type, battery condition and tire pressure.

<sup>4)</sup> With rear caster suspension.

<sup>5)</sup> Does not include seat or battery weights.

<sup>6)</sup> Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
- <sup>7)</sup> 200 lbs. weight capacity.

<sup>8)</sup> Solid and pneumatic tires are available; solid tires may incur an upcharge if not included in package.

<sup>9)</sup> Battery weight may vary based on manufacturer.

<sup>10)</sup> 40A, PG 6-key VR2 required with Power Elevating Seat Option.

<sup>11)</sup> Not available with TRU-Balance<sup>®</sup>

<sup>12)</sup> Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.



facebook.com/PrideMobilityProductsCorp

youtube.com/PrideMobility

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

**Pride Mobility Products Corporation<sup>®</sup>**  
(US) 800-800-8586 • (Canada) 888-570-1113 • [pridemobility.com](http://pridemobility.com)

**Pride<sup>®</sup>**  
Mobility Products Corp.





The Jazzy® Power Chair is engineered and factory-tested to provide you with consistent high-performance operation you can count on. The Jazzy's high-quality construction and proven durability means you can depend on it for a long time.

With several high-performing Jazzy Power Chair models available and many accessories and options to choose from, there's a Jazzy to meet your unique needs and preferences.

Options

- Cup holder
- Oxygen tank holder
- Weather cover
- Cane/crutch holder
- Rear basket
- Walker holder



Jazzy® Elite ES Series

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 300 lbs.  
**Range per Charge:** Up to 15 miles (ES)  
Up to 12 miles (ES-1)  
Up to 11 miles (ES Portable)



Jazzy Red



Jazzy Blue



Jazzy® Select Elite

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 300 lbs.  
**Range per Charge:** Up to 15 miles



Jazzy Red



Jazzy Blue



Jazzy® Select 6

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 300 lbs.  
**Range per Charge:** Up to 15 miles



Jazzy Red



Jazzy Blue



Jazzy® Elite 14

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 300 lbs.  
**Range per Charge:** Up to 9.8 miles



Jazzy Red



Jazzy Blue



Jazzy® 600 ES

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 300 lbs.  
**Range per Charge:** Up to 16.65 miles



Jazzy Red



Jazzy® 1113 ATS

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 300 lbs.  
**Range per Charge:** Up to 10.9 miles



Candy Apple Red



Viper Blue



Onyx Black



Champagne



Jazzy® Sport 2

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 300 lbs.  
**Range per Charge:** Up to 15 miles



Jazzy Red



Jazzy Blue



Jazzy® Elite HD

**Maximum Speed:** Up to 4 mph  
**Weight Capacity:** 450 lbs. (400 lbs. on 18" deep seat configurations)  
**Range per Charge:** Up to 9.8 miles



Jazzy Red



Jazzy Blue



Jazzy® 614 HD

**Maximum Speed:** Up to 4.3 mph  
**Weight Capacity:** 450 lbs.  
**Range per Charge:** Up to 11 miles



Candy Apple Red



Viper Blue



Jazzy® 1450

**Maximum Speed:** Up to 4.5 mph  
**Weight Capacity:** 600 lbs.  
**Range per Charge:** Up to 12.14 miles



Candy Apple Red



Viper Blue