

# Pride® Lifts & Ramps

# PRIDE TOUGH

The Milford Person Lift is a highly adaptable and versatile wheelchair to vehicle transfer solution renowned to be quick and easy to operate for a wide range of people. With a lifting capacity of 330 lbs., the Milford Person Lift can be used anywhere from transferring to your vehicle, in the home and office, or while traveling.

## FEATURES

- Works in a wide variety of vehicles
- Can be fitted to the driver's side, passenger's side, and front, middle or back seats
- 330 lbs. weight capacity
- Sling sizes from extra small to extra large with amputee slings also available
- Lifts quickly and easily with the touch of a button
- Stress-free application for the patient and caregiver
- Emergency switch and optional manual operation



## Milford Person Lift

BASE WEIGHT <sup>3</sup>	18 lbs.
MAXIMUM SCOOTER/PC OVERALL WEIGHT	330 lbs. max.
MAXIMUM PRODUCT WIDTH	N/A
OVERALL PLATFORM LENGTH	N/A
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	22.5"
MOTOR	12 Volt DC
WARRANTY	1-Year Transferable
HITCH CLASS <sup>4</sup>	N/A
OPTIONS	Sling Size & Type, Wall Bracket, Smartbase
VEHICLE PLATFORM	Full-size van, Minivan, Sedan, SUV and Crossover
LIFTS CURRENT PRIDE® PRODUCTS <sup>1,2</sup>	N/A

1. Lifting weight does not include accessories or certain upgraded seating. Please contact Customer Service for greater detail.
2. Vehicle dependent.
3. Does not include any accessories or hardware.
4. Hitch weight carrying capacity varies by vehicle. Call for details.

Truly adjustable transfer solutions, quick & easy to operate



**Milford SmartBase**

Smartbase option available for mobile application

**OPTION**

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —

www.prideliftsandramps.com • www.pridemobility.com  
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)