

Golden Power Wheelchairs

HIGHLY DEPENDABLE WITH MAXIMUM MANEUVERABILITY, SUPERIOR STABILITY AND GREATEST FUNCTIONALITY



Visit Your Local Authorized Dealer:

chnologie



THE MOST STABLE BASE IN ITS CLASS!

YOU CAN'T STAND LIKE THIS ON ANY OTHER POWER WHEELCHAIR BASE

The Golden Compass[™] line of power wheelchairs provide superior forward stability making them the safest power wheelchairs to mount or dismount. The Compass line has fewer moving parts than a power base with springs for the best forward stability and no tipping!

THE INNOVATIVE SEAT DESIGN MAKES OUR POWER WHEELCHAIRS THE EASIEST TO TRANSPORT!

Golden's unique Captain's Seat design found on the Compass Sport, Compass HD and Compass TRO features the patented Fold-Flat Seat Bracket standard. This unique bracket allows the backrest to fold completely flat against the seat, requiring a lower clearance height to load the unit into a vehicle. The Compass HD and TRO also feature the ergonomic Slide Seat Bracket so the seat can be slid off or onto the seat frame for easy assembly or disassembly.

Experience Constant Contact.

Inclines and thresholds are no match for the Golden Compass Sport's true center wheel drive design with the patented springless articulating chassis!

This unique suspension system keeps all six wheels on the ground simultaneously, eliminating the rocking sensation when traversing inclines, giving all Compass[™] models unparalleled curb and threshold climbing ability. The Compass's suspension system pushes the center drive wheel down when the front and rear casters are raised and pushes the front and rear casters down when the drive wheels are raised. No other power base design can match the constant wheel to ground contact of any Golden Compass power base!



Scan QR code to view the Compass product video.

Power Elevating Seat Option:

- Available on GP620
- Seat extends up to 10"
- Can be driven at full speed with seat raised up to 3"
- Can be driven at half speed with seat raised 3"-5"
- Automatic lockout prevents unit from moving when seat is between 5"-10"

Still Love the Thrill of the Chase!

Test drive the industry's most comprehensive line of consumer power wheelchairs today!

- Model GP160 Retail Cash
- **Power Chair!**
- Wt. Cap.: 300 lbs.
- Turning Radius: 26"
- Batteries: 2-18Ah
- Single Post Seat Mount

- Model GP162 LiteRider **ptc** • Retail Cash Power Chair!
 - Wt. Cap.: 300 lbs.
 - Turning Radius: 26"
 - Batteries: 2-18AH
 - Single Post Seat Mount
 - DISASSEMBLES FOR TRANSPORT

- PORT Model GP205 Medicare
 - Codes K0823 & K0822

GOLDEN

- Wt. Cap.: 300 lbs
- Turning Radius: 25"
 - Batteries: 2-U1
 - Single Post Seat Mount



- Model GP207 Medicare Codes: K0816 and K0815
- Wt. Cap.:
- 300 lbs.
- Turning Radius: 25"
- Batteries:
- 2 U-1 • Single Post
- Seat Mount
- tolden golden 1PASS''' • Model GP605 Model GP620 Medicare Codes K0823 & K0822 K0825 & K0824 • Wt. Cap.: 300 lbs. • Wt. Cap.: 450 lbs. • Turning Radius: 19.5" • Turning Radius: 20" Batteries NF22 • Four Post Seat Mount with Slide Seat Bracket Model GP615 • Turning Radius: 20" Medicare Codes • Batteries NF22
- Medicare Codes

golden

- Batteries U1
- Single Post Seat Mount

- K0836 & K0835
- Wt. Cap.: 300 lbs.
- Accepts Many Tilt Seating Systems



POWER WHEELCHAIR PRODUCT SPECIFICATIONS

Model Number		GP605 Compass Sport	GP620 Compass HD	GP615 Compass TRO	GP160 LiteRider PTC	GP162 LiteRider PTC	GP205 Alanté Sport	GP207 Alanté
Medicare Code		K0823/ K0822	K0825/ K0824	K0836/ K0835	Not Coded	Not Coded	K0823/ K0822	K0816/ K0815
Weight Capacity		300 lbs.	450 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.
Drive Wheels		Center	Center	Center	Rear	Rear	Front	Front
Maximum Speed		4 mph	4 mph	4 mph	3.5 mph	3.5 mph	3.5 mph	3.5 mph
Operating Range ¹		15 miles	12 Miles	12 Miles	9.5 miles	9.5 miles	16.5 miles	16.5 miles
Ground Clearance ³		3.25"	3.5"	3.5"	2.5″	2.5″	1.75″	1.75″
Turning Radius		19.5"	20"	20"	26″	26″	25″	25″
Type Batteries Colors		2-U1	2-NF22	2-NF22	2-18 AH	2-18 AH	2-U1	2-U1
		Red, Blue	Red, Blue	Red, Blue	Red	Red	Red	Red
Freewheel Mode		Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electronic Speed Control		Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electro- Mechanical Brakes		Yes	Yes	Yes	Yes	Yes	Yes	Yes
Charger		Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board
Controller		Dynamic Shark 60 Amp Remote	Dynamic Shark 75 Amp Remote	Dynamic Shark 75 Amp Remote	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote
Length		39"	40"	40"	33″	33″	39.5″	39.5″
	Width	24"	25"	25"	22.5″	22.5″	22.25″	22.25″
WARRANTY	Frame	Lifetime	5 Years	5 Years	5 years	5 years	5 years	5 years
	Drive Train	13 months	13 months	13 months	13 months	13 months	13 months	13 months
WA	Electronics ²	13 months	13 months	13 months	13 months	13 months	13 months	13 months

Operating range will vary due to rider weight, drive surface and terrain.
Electronics warranty excludes batteries.

3. Measured at the lowest point mid-frame.

Note: Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.