



GoldenLift Chairs

AMERICA'S #1 SELLING BRAND

FALL 2014 Edition



Visit Your Authorized Golden Dealer:

Our American made lift chairs come in every shape and size for the perfect fit. Our premium fabrics match nearly every décor, and our patented technologies give you the most comfortable lift and recline experience.



MAXICOMFORT® SERIES

(maximum comfort)

MaxiComfort® chairs are the most advanced lift and recline chairs anywhere in the world! In MaxiComfort, your body is transported into a blissful, relaxed state. Choose from optimal comfort positions that promote health, wellness, and relaxation for both your body and your mind.



MAXI COMFORT®
Ultimate Recline Technology

GOLDEN AutoDrive™
CUSTOMIZED COMFORT

Travel effortlessly and smoothly through your favorite positions. Using pre-programmed and programmable buttons, you can lift, recline, watch TV, sleep, or enjoy the wellness benefits of the zero-gravity and Trendelenburg positions, all at the touch of the AutoDrive™ buttons.



When you sit in the Cloud, you feel cushion and support from every angle. Get surrounded with first class comfort in hybrid sizing lift chairs! The Cloud SME comes in small/medium for people 5'1" to 5'6" and the Cloud MLA fits people 5'7" to 6'4".

CLOUD FABRICS



*upcharge applies



MaxiComfort Cloud
PR-510-MLA shown in Shiraz



MaxiComfort Cloud SME
PR-510-SME shown in Calypso



Quick Ship Fabrics for PR-510-SME
Hazelnut, Sterling, Copper

MaxiComfort Positioning



= pocketed coil spring seat std.

MAXICOMFORTER FABRICS



Quick Ship Fabrics for PR-505S, PR-505M & PR-505L: Palomino, Admiral, Cabernet

Quick Ship Fabrics for PR-505JPT: Pearl



MaxiComforter Jr. Petite
PR-505JP shown in Pearl



MaxiComforter Small
PR-505S shown in Cabernet



MaxiComforter Medium
PR-505M shown in Palomino



MaxiComforter Large
PR-505L shown in Admiral



Additional Quick Ship Fabrics for PR-505M: Pearl, Brisa® Buckskin & Coffee Bean (upcharge applies).



Cirrus
PR-508M shown in Palomino

CIRRUS FABRICS



*upcharge applies



Quick Ship Fabrics for PR-508M: Palomino, Admiral, Pearl, and Brisa® Buckskin (upcharge applies).

RELAXER FABRICS



*upcharge applies



Relaxer
PR-756MC shown in Copper



Quick Ship Fabrics for PR-756-MC & PR-756L: Hazelnut, Sterling & Copper.

Additional Quick Ship Fabrics for PR-756-MC: Brisa® Buckskin & Coffee Bean (upcharge applies).



Relaxer
PR-756L shown in Sterling



COMFORTER SERIES

(a perfect fit)

The Comforter series fits people who are 5'0" to 6'8" tall with weight capacities from 300 to 375 lbs. These chairs give the perfect sizes for the perfect fit.



Quick Ship Fabrics for PR-501S, PR-501M and PR-501L: Palomino, Admiral, & Cabernet.



Jr. Petite

PR-501JP shown in Admiral
For heights 5'0" and under



Quick Ship Fabrics for
PR-501JPT: Pearl.



Small

PR-501S shown in Pearl
For heights 5'0" to 5'3"



Medium

PR-501M shown in Palomino
For heights 5'4" to 5'10"



Quick Ship Fabrics for
PR-501M: Pearl.



Large

PR-501L shown in Evergreen
For heights 5'11" to 6'2"



Tall

PR-501T shown in Cabernet
For heights 6'3" to 6'8"




Quick Ship Fabrics
for PR-501T: Palomino.



COMFORTER WIDES

The Comforter Wides are perfect for those who need more space between the armrests. Constructed using the highest quality components for reliable, steady lifting and reclining, these chairs offer weight capacities from 375 to 700 lbs.

COMFORTER WIDE FABRICS



Quick Ship Fabrics for 501M26D & PR-502:
Palomino

PROUDLY



MADE IN AMERICA

with Domestic & Imported Materials



Small 23

PR-501S-23 shown in Cabernet
For heights 5'0" to 5'3"



Medium 26

PR-501M-26 shown in Pearl
For heights 5'4" to 5'10"
(Dual Motor)



Large 26

PR-501L-26 shown in Evergreen
For heights 5'11" to 6'2"
(Dual Motor)



Tall 28

PR-501T-28 shown in Admiral
For heights 6'3" to 6'8"
(Dual Motor)



Super 33

PR-502 shown in Palomino
(Triple Motor)

VALUE SERIES

(cost-effective)

This series combines ultra-soft lumbar pillow support with our unique "grip-rite" arm design for easy transfer.



Monarch Md.
PR-355M shown in Autumn



Capri
PR-200 shown in Cornflower



Monarch Lg.
PR-355L shown in Spring

VALUE FABRICS



Quick Ship Fabrics for PR-200, PR-355M & PR-355L:
Autumn, Cornflower, Rosewood and Chestnut

SIGNATURE SERIES

(multi-functional)

Want a little something extra from your lift chair?
This selection has it all, and... even a little more.



Burlington
PR-722 shown in Calypso
Sink into the comfort of pocketed coil spring seats combined with our plush three-pillow waterfall back and let the stress of the day melt away.



Regal
PR-751 shown in Sterling
Ease into this luxurious 3-pillow lift chair with convenient arm storage and foldable tray.



Savannah
PR-930 shown in Copper
Thanks to our unique frame design, the Savannah can be placed only 3" from the wall and still recline to the full napper position!

SIGNATURE FABRICS



Quick Ship Fabrics for PR-751TY, PR-722, PR-724 & PR-451: Hazelnut
Quick Ship Fabrics for PR-930: Copper
Quick Ship Fabrics for PR-751TY: Sterling & Copper



Shiatsu
PR-724 shown in Sterling
Relax and enjoy the rotating massage mechanism of the Shiatsu lift and recline chair.

TRADITIONAL SERIES

(universal designs)

The Traditional series is all about looks. Choose from colonial or overstuffed biscuit back classics.



Oxford
PR-710 shown in Silt
The classic "club chair" design with a lift! Perfect as a second lift chair in your home, ideal for a bedroom or den, the Oxford features a pocketed coil spring seat and designer nail heads in front of the armrests.

OXFORD FABRICS



W = pocketed coil spring seat std.



Cambridge
PR-401SME shown in Sterling



Cambridge
PR-401MLA shown in Hazelnut



Williamsburg
PR-747 shown in Shiraz

Quick Ship Fabrics for PR-401-SME & 401-MLA:
Hazelnut, Sterling, Copper

CAMBRIDGE & WILLIAMSBURG FABRICS



Transfer Chair
PR-451 shown in Hazelnut
Unique flip-up arm design allows for easy transfer into and out of a power chair.



Golden Technologies is the world's leading manufacturer of lift and recline chairs. Founded in 1985, Golden has become America's most trusted and reliable brand. Our lift chairs are hand assembled in the United States and built using only the highest quality materials.

PRODUCT SPECIFICATIONS

		Model	Weight Capacity	Overall Width	Overall Height	Seat to Top of Back	Width Between Arms	Seat Depth	Floor to Top of Seat	Extension - While Reclined from Back of Seat to End of Footrest	Floor to Top of Footrest - Reclined	Heat & Massage Option Available	Footrest Extension Option Available	Straight Lift Option Available*
Traditional	Cambridge	PR-401-SME	375 lbs.	33"	41.5"	27"	20"	19"	19"	41.5"	22"	Yes	Yes	Yes
	Cambridge	PR-401-MLA	375 lbs.	35"	44.5"	31"	22"	22"	22"	44"	23"	Yes	Yes	Yes
	Williamsburg	PR-747	375 lbs.	32.5"	44"	28"	21"	21.5"	19.5"	42"	23"	Yes	Yes	Yes
	Oxford	PR-710	375 lbs.	38.5"	41"	27"	22"	19.5"	20.5"	43"	22.5"	No	No	Yes
Signature	Regal	PR-751TY	375 lbs.	36"	44"	28"	22"	21.5"	20"	42"	24"	Yes	No	Yes
	Burlington	PR-722	375 lbs.	32.5"	43"	28"	21"	20.5"	20.5"	42.5"	23"	Yes	Yes	Yes
	Shiatsu	PR-724	375 lbs.	32.5"	44"	29"	21"	21"	20"	42"	22.5"	No	No	Yes
	Transfer	PR-451	375 lbs.	33"	42.5"	29"	20"	20"	19"	40.5"	23"	Yes	Yes	Std.
	Savannah	PR-930	350 lbs.	33"	39"	24.5"	17"	18.5"	19"	39.5"	19.5"	No	No	No
MaxiComfort®	Cloud SME	PR-510-SME	375 lbs.	33.5"	41"	29"	18.5"	19.5"	19"	43"	23"	Yes	No	No
	Cloud MLA	PR-510-MLA	375 lbs.	38"	45"	29"	20"	21.5"	20"	44"	24.5"	Yes	No	No
	Relaxer	PR-756MC	375 lbs.	36.5"	44"	29"	20"	21"	19"	43"	23"	Yes	Yes	No
	Relaxer	PR-756L	375 lbs.	39"	46"	30"	22"	22"	20.5"	44"	24"	Yes	Yes	No
	Cirrus	PR-508M	375 lbs.	32"	43"	29"	20"	20"	20.5"	42"	23"	Yes	Yes	No
	Comforter JP	PR-505JP	300 lbs.	30.5"	39"	22.5"	19"	18.5"	18"	38"	22"	Yes	Yes	No
	Comforter Sm.	PR-505S	300 lbs.	31.5"	40"	25"	20"	19.5"	19"	39"	23"	Yes	Yes	No
	Comforter Md.	PR-505M	375 lbs.	32.5"	42"	27"	21"	21"	19.5"	42"	23"	Yes	Yes	No
	Comforter Lg.	PR-505L	375 lbs.	34.5"	44"	28.5"	23"	22"	20.5"	43"	23"	Yes	Yes	No
Comforter	Junior Petite	PR-501JP	300 lbs.	30.5"	39"	25"	19"	18.5"	18"	38"	22"	Yes	No	Yes
	Small	PR-501S	300 lbs.	31.5"	40"	25"	20"	19.5"	19"	39"	23"	Yes	Yes	Yes
	Medium	PR-501M	375 lbs.	32.5"	42"	27"	21"	21"	19.5"	42"	23"	Yes	Yes	Yes
	Large	PR-501L	375 lbs.	34.5"	44"	28.5"	23"	22"	20.5"	43"	23"	Yes	Yes	Yes
	Tall	PR-501T	375 lbs.	34.5"	47"	30"	23"	24"	22.5"	45"	23"	Yes	Yes	Yes
	Small X-wide	PR-501S-23	375 lbs.	34.5"	40"	25"	23"	19.5"	18"	39"	23"	Yes	Yes	Yes*
	Medium-26 Dbl.	PR-501M-26D	500 lbs.	38.5"	42"	27"	27.5"	21.5"	19.5"	42"	23"	Yes	Yes	Yes*
	Large-26 Dbl.	PR-501L-26D	500 lbs.	38.5"	44"	28.5"	27.5"	22.5"	20.5"	43"	23"	Yes	Yes	Yes*
	Tall-28 Dbl.	PR-501T-28D	500 lbs.	41"	48"	30"	29.5"	24"	22.5"	45"	23"	Yes	Yes	Yes*
Value	Super 33	PR-502	700 lbs.	45"	47"	29"	33"	22"	21"	45"	21"	Yes	No	Yes*
	Capri	PR-200	325 lbs.	31"	41"	27"	20"	19"	18"	40"	23"	No	No	Yes
	Monarch Md.	PR-355M	375 lbs.	32"	42"	28"	21"	21.5"	18.5"	43"	23"	Yes	No	Yes
	Monarch Lg.	PR-355L	375 lbs.	34"	44"	30"	23"	22"	18.5"	43"	22"	Yes	No	Yes

*Must be set at factory.

