


MAXIMA

ALSO AVAILABLE IN 3-WHEEL.

COLOR CHOICES



COLOR CHOICES



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

SPECIFICATIONS

	Front Tires	Rear Tires	Maximum Speed <sup>2</sup>	Ground Clearance <sup>4,12,17</sup>	3-W Turning Radius <sup>4</sup>	4-W Turning Radius <sup>4</sup>	3-W Length/Width <sup>4,12</sup>	4-W Length/Width <sup>4,12</sup>	Seat-to-ground Height Range	(3/4-wheel) Weight <sup>2,3,14</sup>	Seat Type (Standard)	Seat-to-deck Height Range	Battery Weight <sup>11</sup>	Battery Requirements <sup>7</sup>	Per Charge Range (Max.) <sup>3,8</sup>	Battery Charger	Maximum Weight Capacity	Weight of Heaviest Piece <sup>4</sup>	Warranty
Go-Go® ES 2 K0800 - GROUP 1 STANDARD	2" x 7.5" SOLID	3" x 7.5" SOLID	UP TO 4 MPH	1.5" at mid-frame	31"	N/A	37"/19.5"	N/A	19"-20"	58.4 LBS.	COMPACT, FOLDABLE (17" X 15")	16.25"-17.25"	6 LBS. EA.	(2)10AH	UP TO 6 MILES	OFF-BOARD, 1.5A	250 LBS.	42 LBS. (SCOOTER BASE)	3-YEAR LIMITED ON FRAME; 1-YEAR LIMITED ON ELECTRONICS; 1-YEAR LIMITED ON DRIVETRAIN
Go-Go® Ultra X	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4 MPH	2"	33"	44"	37"/19.5"	39.5"/19.5"	18.13"-20.13"	75 LBS./ 78.5 LBS.	COMPACT, FOLDABLE (17" X 17")	15.75"-17.25"	9 LBS. EA. <sup>8</sup>	(2)12AH	UP TO 8 MILES	OFF-BOARD, 2A	260 LBS.	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR LIMITED ON ELECTRONICS; 1-YEAR LIMITED ON DRIVETRAIN
Go-Go Elite Traveller® K0800 - GROUP 1 STANDARD	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4 MPH	2"	33"	44"	37"/19.5"	39.5"/19.5"	18.13"-20.13"	75 LBS./ 78.5 LBS.	COMPACT, FOLDABLE (17" X 17")	15.75"-17.25"	9 LBS. EA. <sup>7</sup> / 14 LBS. EA. <sup>8</sup>	(2)12AH OR (2)18AH	3W: UP TO 7.2 MILES <sup>15</sup> 4W: UP TO 7.7 MILES <sup>16</sup>	OFF-BOARD, 2A	300 LBS.	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR LIMITED ON ELECTRONICS; 1-YEAR LIMITED ON DRIVETRAIN
Go-Go Elite Traveller® Plus K0806 - GROUP 2 <sup>11</sup> STANDARD	2" x 8" SOLID	3" x 9" SOLID	UP TO 4.5 MPH	3"	34"	45.75"	40"/21"	40.25"/21"	20.5"-22"	84 LBS./ 86.5 LBS.	COMPACT, FOLDABLE (17" X 17")	16.25"-17.75"	14 LBS. EA. <sup>8</sup>	(2)18AH	3W: UP TO 12.4 MILES 4W: UP TO 13.3 MILES	OFF-BOARD, 2A	300 LBS.	32.5 LBS. (3-WHEEL) (FRONT SECTION) 35 LBS. (4-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR LIMITED ON ELECTRONICS; 1-YEAR LIMITED ON DRIVETRAIN
Go-Go® Sport K0801 - GROUP 1 HEAVY DUTY	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	3" x 9" SOLID	UP TO 4 MPH	3"	36"	52"	43"/20"	43"/20"	20.5"-22"	83.2 LBS./ 99.2 LBS.	COMPACT, FOLDABLE (17" X 17")	16.25"-17.75"	14 LBS. EA. <sup>8</sup>	(2)18AH	3W: UP TO 9 MILES 4W: UP TO 12 MILES	OFF-BOARD, 2A	325 LBS.	36.2 LBS. (3-WHEEL) (FRONT SECTION) 52.2 LBS. (4-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR LIMITED ON ELECTRONICS; 1-YEAR LIMITED ON DRIVETRAIN
Go-Go® LX with CTS Suspension	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4.4 MPH	1.5" at motor	33"	44.75"	39.75"/20.5"	41"/20.5"	20.75"-22.25"	85.4 LBS./ 87.4 LBS.	COMPACT, FOLDABLE (17" X 17")	16.25"-17.25"	9 LBS. EA. <sup>7</sup> / 14 LBS. EA. <sup>8</sup>	(2)12AH OR (2)18AH	3W: UP TO 9.44 MILES 4W: UP TO 10.5 MILES	OFF-BOARD, 2A	300 LBS.	34 LBS. (3-WHEEL) (FRONT SECTION) 36.6 LBS. (4-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR LIMITED ON ELECTRONICS; 1-YEAR LIMITED ON DRIVETRAIN
Victory® 9 K0800 - GROUP 1 STANDARD	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	2.5" x 9" SOLID	UP TO 5.25 MPH	3W: 2" 4W: 2.25"	40.25"	51"	43.5"/22.25"	45.5"/22.5"	19.25"-20.75"	117.5 LBS./ 131 LBS.	FOLDABLE, MOLDED PLASTIC (18"X17")	17.25"-18.75"	24.5 LBS. EA.	(2) U1	UP TO 13 MILES	OFF-BOARD, 3A	300 LBS.	48.5 LBS. (3-WHEEL) (FRONT SECTION) 59 LBS. (4-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR PRO-RATA ON DRIVETRAIN; 1-YEAR LIMITED ON ELECTRONICS
Victory® 10 K0801 - GROUP 1 HEAVY DUTY	3.5" x 10" SOLID	3.5" x 10" SOLID	UP TO 5.25 MPH	3" at mid-frame	42"	54"	45.25"/22.25"	47"/22.25"	20.5"-22"	131.5 LBS./ 132.5 LBS.	FOLDABLE, MOLDED PLASTIC (18"X17")	17.5"-19"	24.5 LBS. EA./ 32.5 LBS. EA.	(2) U1 OR (2) 40AH	UP TO 15.5 MILES	OFF-BOARD, 3A	400 LBS.	54.5 LBS. (3-WHEEL) (FRONT SECTION) 61 LBS. (4-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR PRO-RATA ON DRIVETRAIN; 1-YEAR LIMITED ON ELECTRONICS
Victory® Sport	3.5" x 10" SOLID	4" x 10" SOLID	UP TO 8 MPH	2.25"	N/A	51.75"	N/A	47"/26"	21.5"-24.25"	150 LBS.	CRS RECLINING H/B W/SLIDERS & HEADREST (18"X18")	17.25"-18.75"	32.5 LBS. EA.	(2) 40AH	UP TO 11 MILES	OFF-BOARD, 5A	350 LBS.	68 LBS. (3-WHEEL) (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR PRO-RATA ON DRIVETRAIN; 1-YEAR LIMITED ON ELECTRONICS
Maxima K0802 - GROUP 1 VERY HEAVY DUTY (3W) K0808 - GROUP 2 <sup>11</sup> VERY HEAVY DUTY (4W)	3" X 10" SOLID	4" X 10" SOLID	UP TO 5.25 MPH	3.25"	49.25"	54"	46"/25.5"	46.9"/22.5"	22.5"-24.5"	186 LBS./ 192 LBS.	DELUXE CONTOURED MEDIUM-BACK (22"X18")	16"-18"	U-1: 24.5 LBS. EA./ NF-22: 38 LBS. EA.	(2) U1 OR (2) NF-22	3W: UP TO 15.4 MILES/NF-22 3W: UP TO 8.7 MILES/U-1 4W: UP TO 16 MILES/ NF-22 4W: UP TO 10.2 MILES/U-1	OFF-BOARD, 5A	500 LBS.	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR PRO-RATA ON DRIVETRAIN; 1-YEAR LIMITED ON ELECTRONICS
Pursuit®	3" X 13" SOLID	13" SOLID	UP TO 7.25 MPH	3.56"	N/A	67.25"	N/A	53.25"/25.25"	24.75"-26.25"	218.5 LBS.	CRS RECLINING H/B W/SLIDERS & HEADREST (19"X17")	17"-18.75"	38 LBS. EA.	(2) NF-22	UP TO 14 MILES	OFF-BOARD, 5A	400 LBS.	95 LBS. (FRONT SECTION)	3-YEAR LIMITED ON FRAME; 1-YEAR PRO-RATA ON DRIVETRAIN; 1-YEAR LIMITED ON ELECTRONICS
Pursuit® XL	3" X 13" PNEUMATIC	13" PNEUMATIC	UP TO 8.5 MPH	5.25"	N/A	82.5"	N/A	56.25"/27.25"	26.25"-27.75"	259.9 LBS.	CRS RECLINING H/B W/SLIDERS & HEADREST (18"X18")	21"-22.5"	60.5 LBS. EA.	(2) GROUP 27 OR (2) GROUP 24	G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES	OFF-BOARD, 8A	400 LBS.	205.5 LBS. (MAIN FRAME)	3-YEAR LIMITED ON FRAME; 1-YEAR PRO-RATA ON DRIVETRAIN; 1-YEAR LIMITED ON ELECTRONICS

1. All shroud panels included.

2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

3. Battery weight not included.

4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

5. Battery pack is 18 lbs. with 12AH batteries.

6. Durable solid plastic.

7. Deep cycle (AGM or gel cell type recommended).

8. 18 AH battery pack weighs 28 lbs.

9. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance.

Test conducted at maximum weight capacity.

10. Up to 14 miles per charge on 18 AH batteries.

11. Battery weight may vary based on manufacturer.

12. Width of base (not including seat).

13. Does not account for anti-tip clearance.

14. Includes standard seat weight.

15. Requires upgrade to medium or high-back contoured seat; additional charges apply.

16. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.

17. At midframe.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.



SCGUIDE-04-07-14



THE ULTIMATE IN STYLE  
& PERFORMANCE®

PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

MOBILITY SCOOTER GUIDE

www.pridemobility.com • www.pridescooters.com

www.pridemobility.com • www.pridescooters.com



COLOR CHOICE



Black



ALSO AVAILABLE IN 3-WHEEL

COLOR CHOICES



Red



Blue



Silver

INTERCHANGEABLE SHROUDS<sup>†</sup>



ALSO AVAILABLE IN 4-WHEEL

COLOR CHOICES



Red



Blue

INTERCHANGEABLE SHROUDS<sup>†</sup>



ALSO AVAILABLE IN 3-WHEEL

COLOR CHOICES



Candy Apple Red



Viper Blue



ALSO AVAILABLE IN 3-WHEEL

COLOR CHOICES



Red



Blue



ALSO AVAILABLE IN 3-WHEEL

COLOR CHOICES



Red



Blue

INTERCHANGEABLE SHROUDS<sup>†</sup>



ALSO AVAILABLE IN 4-WHEEL

COLOR CHOICES



Candy Apple Red



Viper Blue



COLOR CHOICES



Candy Apple Red



Viper Blue

ALSO AVAILABLE IN 4-WHEEL

COLOR CHOICES



Red



Blue

INTERCHANGEABLE SHROUDS<sup>†</sup>

