

Safe Independent Mobility



Shown with optional Seat Base Covers

Seat Bases & Covers are available
in both Gray and Tan

Features and Benefits:

- New 24 Month Warranty on Parts and 12 Month Limited Warranty on Labor.
- New 500 lbs total combined weight capacity eliminates the need for both a standard and a heavy-duty model.
- Test certification to meet the combined FMVSS and CMVSS 207 & 210 requirement when the seatbase is in the full-forward position.
- Improved Rack & Pinion design provides quieter, faster horizontal travel.
- New Central Wiring Distribution Channel.
- New repositioned Wiring Harness eliminates the opportunity for accidental damage.
- New Direct Vehicle Chassis Ground System.
- New 12" diameter Swivel Gear provides smoother turning.
- New Swivel Adjustment Screw allows the seat to be adjusted up to 5° towards center of vehicle.
- New Delrin Bearings provide improved stability and quieter operation.
- New Actuator Protective Cover prevents objects from falling into the actuator.
- New Two Position Adjustable Switch Box allows you to position it as close to the seat as possible.
- New application specific part numbers.
- Signature B&D Direct Mounting System: No Drilling Required - On most models.


www.bdindependence.com

Leadership Series Transfer Seat Base

B & D Independence, Inc. Limited Warranty

Transfer Seat Bases manufactured by B&D Independence, Inc. are guaranteed to be free from defects in material and workmanship for a period of 2 years from the date the product is placed in service. B&D Independence will repair or replace, at our option, any factory installed part that is found to be defective. For a period of 1 year, limited labor charges may be covered by this warranty. Warranty claims must be submitted to B&D Independence during the warranty period by the installing dealer. Some limitations apply.

Authorized Dealer



This Authorized Dealer meets our high standards for sales and service and is qualified to install and service B&D Independence Power Transfer Seat Bases. Authorized Dealers are committed to providing the highest quality service and will honor all manufacturers' warranties.

B&D
Independence, Inc.

Safe Independent Mobility

www.bdindependence.com

B&D
Independence, Inc.

***B&D reliability
has never
looked this good!***



www.bdindependence.com

We were founded in 1981 on the premise that quality and service would never be compromised. These basic principles have served us well, and today B&D Independence, Inc. is the recognized industry leader in Transfer Seat Bases.

Mission Statement:

We are committed to:

- Designing and manufacturing the industry's best transfer seat base at a fair value to our dealers and those who use our products.
- Fully testing each new design to insure it meets or exceeds the FMVSS and CMVSS 207 and 210 combined requirements. (See Statement Below)
- Listening to our customers needs and constantly looking for ways to improve our products.
- Shipping 100% of our shipments within 24 hours.
- Providing competitive wages and benefits for our employees and recognizing their contribution in helping to achieve our vision.
- Treating everyone we come in contact with each day with respect and dignity.

Our Vision


We will be recognized throughout the world as the mobility industry leader for quality, safety and customer service.



Transfer Seat Bases manufactured by B & D Independence, Inc. meet the combined Federal Motor Vehicle Safety Standards (FMVSS 207 & 210) and Canadian Motor Vehicle Safety Standards (CMVSS 207 & 210) when installed according to the B&D Independence installation instructions and driven in the full forward position.

Leadership Series Transfer Seat Base

Model L51



6 way seat base for easy transfer



Models available for most mini & full sized vans
See dealer for details



New from B&D

"Peace of Mind" Mechanical Release

We are pleased to announce that B&D Independence transfer seat bases for lowered floor minivans are available with a mechanical release that we call "Peace of Mind". This release allows you to move the seat base away from the steering wheel in the event of a power failure.

The "Peace of Mind" Mechanical Release is mounted below the switch box controls, and is conveniently located for easy access for the driver, passenger or caregiver. It is designed to prevent accidentally engaging the mechanical release, while at the same time being accessible and easy to operate.

Now you can have "Peace of Mind" and America's most popular transfer seat base!

B&D reliability has never looked this good!

Powder Coat Options

Approximately 85% of vans have gray and/or charcoal interior. For these vehicles, we offer a medium gray powder coat for a long lasting, quality finish. Of the remaining 15% of vans, most have a tan or champagne interior. We are pleased to announce that, with the addition of an enhanced powder system, we now have the ability to powder coat our seat bases tan. Contact your mobility dealer for all the details.

Seat Base Covers

The appearance of the equipment is second only to your safety. That's why we have developed our new ABS seat base covers that are now available for most lowered floor minivans. The 3 piece design comes with all the hardware needed for your dealer to install and adds the finishing touch to America's most popular transfer seat base. These OEM quality covers are available in gray and tan for both driver and passenger positions. Your mobility dealer can provide you with all the details.